


WASHER PROCESS MONITORING SYSTEM

CENTRAL STERILE PROGRAM



A comprehensive verification and data management solution for automated and ultrasonic cleaning equipment

WASHER PROCESS INDICATOR & HOLDER

- ▲ Sensitive to aspects of the automated and ultrasonic cleaning process and validated against a standard test soil*
- ▲ Mimics the realistic cleaning challenge posed by hard-to-reach instrument areas
- ▲ No blood products and non-toxic residuals
- ▲ Easy to interpret test soil with bright color and high contrast

WASHER PROCESS MANAGER

- ▲ Captures verification test results and data using an iOS app
- ▲ Customizable to meet your facility standard practices
- ▲ Automated cloud-based data management

ONLINE DATA PORTAL

- ▲ Real-time dashboard and on-demand data analysis
- ▲ Powerful insights for troubleshooting and staff training
- ▲ Customizable to meet facility-specific needs

Test results from Washer Process Indicators are digitally captured on the Washer Process Manager App, and stored on the cloud-based customer Data Portal for on-demand data analysis and reporting



Talk to Your Representative About Our **Ecolab® Central Sterile Program** featuring the Washer Process Monitoring System

Part #	Description	Qty/Case
92023641	Washer Process Indicator	200
92023642	Washer Disinfectant Process Indicator Holder	1
53000615	Ultrasonic Washer Process Indicator Holder	1

NOTE: The Ecolab Washer Process Monitoring System is only available as a part of the Ecolab Central Sterile Program

*Soil validated against ISO/TS 15883-5:2005 Standard Annex Q